


iDOL

CERRAHI ALETLER
WWW.idolcerrahialetler.com

Arthroscopy Instruments





Mod. PAA
















Mod. PAB



Mod. PAC














Arthroscopy

All Instruments with Overload Protection

		Hook Scissors 	Hook Punch with Scoop  L = 5,0 mm · B = 1,6 mm	Hook Punch set down  L = 5,0 mm · B = 1,6 mm	Hook Punch Cut angular  L = 5,0 mm · B = 1,6 mm
straight 	2,7	xxx 250-27	xxx 290-27	xxx 350-27	
	3,5	xxx 250-35	xxx 290-35	xxx 350-35	xxx 370-35
30° right 	2,7	xxx 251-27	xxx 291-27	xxx 351-27	
	3,5	xxx 251-35	xxx 291-35	xxx 351-35	xxx 371-35
30° left 	2,7	xxx 253-27	xxx 293-27	xxx 353-27	
	3,5	xxx 253-35	xxx 293-35	xxx 353-35	xxx 373-35
45° right 	2,7	xxx 255-27	xxx 295-27	xxx 355-27	
	3,5	xxx 255-35	xxx 295-35	xxx 355-35	xxx 375-35
45° left 	2,7	xxx 257-27	xxx 297-27	xxx 357-27	
	3,5	xxx 257-35	xxx 297-35	xxx 357-35	xxx 377-35
7° up 	2,7	xxx 259-27	xxx 299-27	xxx 359-27	
	3,5	xxx 259-35	xxx 299-35	xxx 359-35	xxx 379-35
7° down 	2,7	xxx 261-27	xxx 301-27	xxx 361-27	
	3,5	xxx 261-35	xxx 301-35	xxx 361-35	xxx 381-35
15° up 	2,7	xxx 263-27	xxx 303-27	xxx 363-27	
	3,5	xxx 263-35	xxx 303-35	xxx 363-35	xxx 383-35
15° down 	2,7	xxx 265-27	xxx 305-27	xxx 365-27	
	3,5	xxx 265-35	xxx 305-35	xxx 365-35	xxx 385-35














Arthroscopy






All Instruments with Overload Protection



		Hook Punch 15° up milled  L = 5,0 mm · B = 1,6 mm	Wide Jaw Punch large  L = 5,0 mm · B = 4,5 mm	Wide Jaw Punch medium  L = 5,0 mm · B = 3,5 mm	Wide Jaw Punch small  L = 6,0 mm · B = 2,5 mm
	3,5	xxx 390-35	xxx 490-35	xxx 560-35	xxx 630-35
	3,5	xxx 391-35	xxx 491-35	xxx 561-35	xxx 631-35
	3,5	xxx 393-35	xxx 493-35	xxx 563-35	xxx 633-35
	3,5	xxx 395-35	xxx 495-35	xxx 565-35	xxx 635-35
	3,5	xxx 397-35	xxx 497-35	xxx 567-35	xxx 637-35
	3,5	xxx 399-35	xxx 499-35	xxx 569-35	xxx 639-35
	3,5	xxx 401-35	xxx 501-35	xxx 571-35	xxx 641-35
	3,5	xxx 403-35	xxx 503-35	xxx 573-35	xxx 643-35
	3,5	xxx 405-35	xxx 505-35	xxx 575-35	xxx 645-35

Arthroscopy

All Instruments with Overload Protection














		Wide Jaw Punch large 15° up milled  L = 5,0 mm - B = 4,5 mm	Wide Jaw Punch medium 15° up milled  L = 5,0 mm - B = 3,5 mm	Wide Jaw Punch small 15° up milled  L = 6,0 mm - B = 2,5 mm	Wide Jaw Punch small 45° up milled  L = 6,0 mm - B = 2,0 mm
straight 	3,5	xxx 530-35	xxx 600-35	xxx 680-35	xxx 698-35
30° right 	3,5	xxx 531-35	xxx 601-35	xxx 681-35	
30° left 	3,5	xxx 533-35	xxx 603-35	xxx 683-35	
45° right 	3,5	xxx 535-35	xxx 605-35	xxx 685-35	
45° left 	3,5	xxx 537-35	xxx 607-35	xxx 687-35	
7° up 	3,5	xxx 539-35	xxx 609-35	xxx 689-35	
7° down 	3,5	xxx 541-35	xxx 611-35	xxx 691-35	
15° up 	3,5	xxx 543-35	xxx 613-35	xxx 693-35	
15° down 	3,5	xxx 545-35	xxx 615-35	xxx 695-35	

		45° Hook Scissors	45° Punch with Scoop	45° Punch blunt	90° Punch with Scoop	90° Punch blunt
			 L = 4,5 mm - B = 1,6 mm	 L = 4,5 mm - B = 1,6 mm	 L = 4,5 mm - B = 1,6 mm	 L = 4,5 mm - B = 1,6 mm
right angled	3,5	xxx 270-35	xxx 410-35	xxx 420-35	xxx 450-35	xxx 460-35
left angled	3,5	xxx 280-35	xxx 430-35	xxx 440-35	xxx 470-35	xxx 480-35

		90° Wide Jaw Punch medium	90° Wide Jaw Punch small
		 L = 7,0 mm - B = 3,0 mm	 L = 7,0 mm - B = 1,8 mm
right angled	3,5	xxx 740-35	xxx 760-35
left angled	3,5	xxx 750-35	xxx 770-35

Arthroscopy




All Instruments with Overload Protection

		Scissors Punch  L = 5,0 mm - B = 0,8 mm	Scissors Punch 15° up milled  L = 5,0 mm - B = 0,8 mm	Punch Retrograde Cut left 	Punch Retrograde Cut right 
straight 	3,5	xxx 700-35	xxx 720-35	xxx 780-35 L = 4,0 mm - B = 2,5 mm	xxx 790-35 L = 4,0 mm - B = 2,5 mm
30° right 	3,5	xxx 701-35	xxx 721-35		
30° left 	3,5	xxx 703-35	xxx 723-35		
45° right 	3,5	xxx 705-35	xxx 725-35		
45° left 	3,5	xxx 707-35	xxx 727-35		
7° up 	3,5	xxx 709-35	xxx 729-35		
7° down 	3,5	xxx 711-35	xxx 731-35		
15° up 	3,5	xxx 713-35	xxx 733-35		
15° down 	3,5	xxx 715-35	xxx 735-35		

Arthroscopy





All Instruments with Overload Protection

		Grasping Forceps 	Grasping Forceps 1 x 2 Teeth, fine 	Grasping Forceps 1 x 2 Teeth, crass 
	2,7	xxx 800-27		
	3,5	xxx 800-35	xxx 820-35	xxx 840-35

		Grasping and Biopsy Forceps 	Biopsy Spoon Forceps 	Bulldog -Grasping Forceps 1x2 teeth, crass 
	2,7			
	3,5	xxx 860-35	xxx 880-35	
	4,5			xxx 845-35

Arthroscopy

All Instruments with Overload Protection

	Suture Grasping Forceps	Suture-Pierce 45° sharp up milled	Suture-Pierce 10° sharp up milled	Suture-Pierce sharp straight
				
	xxx 890-35	xxx 893-35	xxx 895-35	xxx 900-35



PP 180-01
Caroli Tendon Grasping Forceps



PP 180-00

Tendon Grasping Forceps with ergonomic handle
straight
with ratchet
overall length 190 mm



PP 180-02

Tendon Grasping Forceps
right
with ratchet



PP 180-03

Tendon Grasping Forceps
left
with ratchet



With ergonomic handle

PP 170-00

Mini Grasping Forceps 1x2 teeth with ratchet
Working length: 130 mm



PP 190-00

Meniscus grasping forceps, strong, 4.5 mm x 15 mm jaws length
Working length: 120 mm

New Arthroscopy Collection

Advantages of the new series: High cutting power with light pressure, Optimal precise cutting , at high load capacity.



Mod. PAD



Mod. PAG



Wide jaw Punch

- x 320-01 straight 3,2 mm
- x 320-02 straight 3,2 mm , 15° up milled
- x 320-03 sheath upswept 3,2 mm
- x 320-04 upbiter right, 3,2 mm
- x 320-05 upbiter left , 3,2 mm



Grasping Forceps cupped
x 350-01 straight 3,4mm



Alligator Forceps
x 350-02 straight 3,0 mm



Mod. PAD



Hook Punch with scoop
 Retrograd Punch

- x 160-01 straight 1,7 mm
- x 200-01 straight 2,0 mm
- x 160-02 sheath upbiter 1,6 mm
- x 160-03 straight 1,6 mm 15° up milled
- x 160-04 right curved 1,6 mm
- x 160-05 left curved 1,6 mm



Scissor Punch

- x 140-01 straight 1,4mm
- x 140-02 upbiter 1,4mm
- x 140-03 right curved 1,4mm
- x 140-04 left curved 1,4mm



Hook Punch blunt

- x 160-06 straight 1,6 mm
- x 160-07 sheath upbiter 1,6 mm
- x 160-08 straight 1,6mm,15° up milled
- x 160-09 right curved 1,6mm
- x 160-10 left curved 1,6mm



Mod. PAG



- x 240-01 right cutting
- x 240-02 left cutting



Punch oval

- x 440-01 straight 4,4 mm
- x 440-02 upbiter 4,4 mm
- x 440-03 right curved 4,4 mm
- x 440-04 left curved 4,4 mm



Punch 90°

- x 340-01 right cutting
- x 340-02 left cutting