

The background of the page features several surgical instruments. On the left, a long, thin pair of forceps extends from the top to the bottom. On the right, there is a cluster of various surgical tools, including curved forceps and a grasper, arranged in a fan-like pattern.

VATS
(Video-Assisted Thoracic Surgery)
Brochure

Duval Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 100-01

curved left,
triangular jaw,
25 mm

24 cm

Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 100-02

straight, oblong
ring jaw,
11 mm, slim

24 cm

Code	Description	Working Length
------	-------------	----------------

VT 101-01	curved left, oval ring jaw, 12 mm slim	24 cm
-----------	--	-------

Foerster Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 101-02	curved left, oval ring jaw, 20 mm	24 cm
-----------	---	-------

Foerster Grasping Clamp



Duval Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 101-03

curved left,
triangular jaw,
15 mm

24 cm

Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 101-04

curved left,
oblong ring jaw,
11 mm, slim

24 cm

Allis Grasping Clamp



Code	Description	Working Length
------	-------------	----------------

VT 101-08

curved left,
5x6 toothed jaw,
6 mm, slim

24 cm

Code	Description	Working Length
------	-------------	----------------

VT 100-09	curved left, tapered jaw, slim	24 cm
-----------	--------------------------------------	-------

**Viper Grasping and
Dissecting Clamp**



Code	Description	Working Length
------	-------------	----------------

VT 100-10	curved, DeBakey jaw, slim	24 cm
-----------	---------------------------------	-------

Mixer Dissecting Clamp



Dennis Dissecting Clamp



Code	Description	Working Length
VT 101-05	curved left, slim	24 cm

Dissecting Clamp



90°

Code	Description	Working Length
VT 101-07	DeBakey/Cooley jaw, angled 90°	24 cm

Code	Description	Working Length
------	-------------	----------------

VT 101-09	curved, oval ring jaw, 6 mm, slim	28 cm
-----------	---	-------

Biopsy Forceps




Code	Description	Working Length
------	-------------	----------------

VT 101-10	curved, oval cupped jaw, 7 mm, slim	28 cm
-----------	---	-------

Biopsy Forceps



VT 100-05



VT 100-03



Code	Description	Working Length
VT 100-05	straight, 1 x 2 DeBakey jaw, 100 mm	27 cm
VT 100-03	curved left, 1 x 2 DeBakey jaw, 50 mm	21 cm
VT 100-04	curved left, 1 x 2 DeBakey jaw, 75 mm	23,5 cm
VT 101-06	curved left, 2 x 3 DeBakey jaw, 100 mm	26 cm

VT 100-04



VT 101-06



Harken Clamp


VT 100-06



Code	Description	Working Length
VT 100-06	curved, 1 x 2 DeBakey jaw, 95 mm	24 cm
VT 100-08	curved, 1 x 2 DeBakey jaw, angled shank	21 cm
VT 100-07	curved, 1 x 2 DeBakey jaw, 140 mm	27 cm

Lambert-Kay Clamp

VT 100-08



Lambert-Kay Harken Clamp

VT 100-07



Needle Holder



Code	Description	Working Length
VT 103-01	straight, slim	21,5 cm

Needle Holder



Code	Description	Working Length
VT 103-02	curved, slim	21,5 cm

Metzenbaum Scissors



Code	Description	Working Length
VT 102-01	curved, blunt jaw	19,5 cm

Code	Description	Working Length
VT 104-01	curved shaft	27 cm

Tube Passer Forceps

